




Medical Trolleys







N6 Series-Medical Trolley







N6 series trolley define a new standard of organization and function with specialised accessories systems for the healthcare sector. This series offers stable, smooth mobility and enhanced convenience with changeable accessories system. As an innovation product, the design is based on professional marketing research, customer feedback and healthcare requirements. For releasing patients and nurses stress, the trolley body provides bright color, which changes the traditional depressing impression of hospital and clinics.

Description:



- 1.Surface: **ABS molded table mounting stainless steel plate**, features abrasion resistance and corrosion resistance; three-sided **stainless steel fence** to avoid any falling of materials.
- 2.Pillar: **Aluminum alloy pillar** with powder coated.
- 3.Board: High-quality **aluminum composite panel for side board and back board**.
- 4.Drawer: All drawers are detachable design, easy cleaning and replace. **Zinc coated steel drawer**.
- 5.Drawer handle: ABS
- 6.Spaces: Lager and classified storage with **PC transparent dividers Or ABS Medicine Box**.
- 7.Castor: Noise-free caster with brake, of which diameter is 100mm.
- 8.Round bumper: on the bottom of trolley body.
- 9.Loading: Drawer loading: 50KG; Overall loading: 200KG.

Product name	Pictures	Model No.	Dimensions (mm)		Standard Accessories:
Emergency Trolley		JDEQJ264A	560x475x1000		1. CRS Defibrillator Shelf; 2. Extension Table; 3. 'L'Shape I.V. Pole; 4. CPR Board; 5. Waste Bin; 6. External Power Holder And Plug Board; 7. Oxygen Tank Holder; 8. PC Transparent Divider; 9. SS Wire Basket; 10. Castor With Brake and Directional; 11. Non-slip Mat; 12. Soft Closing Drawer; 13. Seals or Central Locks.
		JDEQJ234A	650x475x1000		
		JDEQJ294A	750x475x1000		
		JDEQJ264B	560x475x970		
		JDEQJ234B	650x475x970		
		JDEQJ294B	750x475x970		
		JDEQJ264D	560x475x970		
		JDEQJ234D	650x475x970		
		JDEQJ294D	750x475x970		

Product name	Pictures	Model No.	Dimensions (mm)		Standard Accessories:
Anesthesia Trolley		JDEMZ264A	560x475x970		1. Extension Table; 2. Waste Bin; 3. PC Transparent Divider; 4. SS Wire Basket; 5. Castor With Brake and Directional; 6. Non-slip Mat; 7. Soft Closing Drawer ; 8. Central Lock.
		JDEMZ234A	650x475x970		
		JDEMZ294A	750x475x970		
		JDEMZ264B	560x475x970		
		JDEMZ234B	650x475x970		
		JDEMZ294B	750x475x970		
		JDEMZ264C	560x475x970		
		JDEMZ234C	650x475x970		
		JDEMZ294C	750x475x970		
Medicine Trolley	 A	JDEFY264A	560x475x900		1. Extension Table; 2. Waste Bin; 3. PC Transparent Divider; 4. SS Wire Basket; 5. Castor With Brake and Directional;; 6. Non-slip Mat; 7. Soft Closing Drawer; 8. Central Lock.
		JDEFY234A	650x475x900		
		JDEFY294A	750x475x900		
	 R	JDEFY264B	560x475x900		
		JDEFY234B	650x475x900		
		JDEFY294B	750x475x900		
		JDEFY294C	750x475x900		

Product name	Pictures	Model No.	Dimensions (mm)		Standard Accessories:
Treatment Trolley		JDEZL264A	560x475x900		1. Extension Table; 2. Waste Bin; 3. PC Transparent Divider; 4. SS Wire Basket; 5. Castor With Brake and Directional; 6. Non-slip Mat; 7. Soft Closing Drawer; 8. Central Lock(1 -drawers trolley not available).
		JDEZL234A	650x475x900		
		JDEZL294A	750x475x900		
		JDEZL264B	560x475x900		
		JDEZL234B	650x475x900		
		JDEZL294B	750x475x900		
		JDEZL264C	560x475x900		
		JDEZL234C	650x475x900		
		JDEZL294C	750x475x900		
		JDEZL264D	560x475x900		
		JDEZL234D	650x475x900		
		JDEZL294D	750x475x900		
Infusion Trolley		JDESE264A	560x475x900		1. Extension Table; 2. Waste Bin; 3. PC Transparent Divider; 4. SS Wire Basket; 5. Castor With Brake and Directional; 6. Non-slip Mat; 7. Needle Box; 8. Soft Closing Drawer;
		JDESE234A	650x475x900		
		JDESE294A	750x475x900		
		JDESE264B	560x475x900		
		JDESE234B	650x475x900		
		JDESE294B	750x475x900		

Product name	Pictures	Model No.	Dimensions (mm)		Standard Accessories:
		JDESE264C	560x475x900		9. Central Lock.
		JDESE234C	650x475x900		
		JDESE294C	750x475x900		
Medical Record Trolley		JDEBL234A	650x475x1070		1.Soft Closing Drawer; 2.Castor With Brake With Brake and Directional; 3.Non-slip Mat.
		JDEBL234B	650x475x1150		
		JDEBL234C	650x475x1010		
Difficult Airway Trolley		JDEQD234A JDEQD234B	800x475x970/1065 445x110x960		1. Extension Table; 2. Central Locks; 3. Castor With Brake With Brake and Directional; 4. Non-slip Mat; 5. Soft Closing Drawer 6.Scopy Box
Housekeeping Trolley	 B	JDEHL264A	1050x475x900		1. Soft Closing Drawer; 2. Castor With Brake With Brake and Directional; 3. Non-slip Mat
	 C	JDEHL264B	1050x475x900		

Product name	Pictures	Model No.	Dimensions (mm)		Standard Accessories:
Instrument Trolley		JDEQJ264A	560x475x970		1.Castor With Brake and Directional.
		JDEQJ234A	650x475x970		
		JDEQJ294A	750x475x970		
		JDEYQ264B	560x475x970		
		JDEYQ234B	650x475x970		
		JDEYQ294B	750x475x970		