

BD BACTEC™ FX40

Blood Culture System

Expanding on the leadership of the BD BACTEC™ FX Blood Culture System.



Improving Workflow Efficiencies...

Simple rapid workflow

- Vial activation
- Enhanced visual indicators
- Instant bottle and instrument status recognition
- Computer touch screen


Optimized Ergonomics

- Modular design provides flexible lab layout

From Specimen Collection to Actionable Results

BD BACTEC Blood Culture Solutions provide proven performance to enhance clinical outcomes and optimize efficiency.

Modular design allows system to grow with your lab demands, additional systems **can easily be added on** as an institution's blood culture volume grows.



Key features of the BD BACTEC™ FX40

LIS connectivity	✓
BD EpiCenter™ connectivity	✓
Blood Volume Monitoring*	✓
Satellite Blood Culturing**	✓
Proven, robust design	✓

* Requires BVM software with EpiCenter

** Requires proper connectivity

Physical Dimensions: Single Instrument

Height	39.1 cm (15.4 in)
Width	67.5 cm (26.6 in)
Depth	58.5 cm (23.0 in)
Clearance (rear, left, right)	0 cm, 7 cm, 38 cm*
Clearance (front)	61.0 cm (24 in)
Weight (empty)	31.8 kg (70.1 lb)

* If tablet is mounted

BD BACTEC™ Plus Aerobic and Anaerobic Lytic media provide ideal recovery of microorganisms.³

